

C11orf48 antibody

Catalog Number: 108637

Product name

C11orf48 antibody

Specificity

Human; other species not tested.

Antibody description

C11orf48 Rabbit Polyclonal antibody. Positive IP detected in HepG2 cells. Positive WB detected in Transfected HEK-293 cells. Observed molecular weight by Western-blot: 28-32 kDa

Preparation

This antibody was obtained by immunization of C11orf48 recombinant protein (Accession Number: NM_001394611). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108637(C11orf48 Antibody) at dilution of



1:700

IP Result of anti-C11orf48 (IP:Catalog No:108637, 3ug; Detection:Catalog No:108637 1:500) with HepG2 cells lysate 1440ug.