

# CXCL12 antibody



Catalog Number: 109637

## Product name

CXCL12 antibody

ELISA, IHC

## Specificity

Human, Mouse; other species not tested.

## Dilutions

Recommended Dilution:

IHC: 1:20-1:200

## Antibody description

CXCL12 Rabbit Polyclonal antibody. Positive IHC detected in human liver cancer tissue, human tonsillitis tissue.

## Preparation

This antibody was obtained by immunization of CXCL12 recombinant protein (Accession Number: NM\_199168). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality


Polyclonal

## Ig Type


Rabbit IgG

## Applications

## Validations



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:109637(CXCL12 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:109637(CXCL12 Antibody) at dilution of 1:50 (under 40x lens)