

ZNF415 Polyclonal Antibody



Catalog Number: 162337

Product name

ZNF415 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ZNF415

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 100-260 of human ZNF415 (NP_060825.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Immunohistochemistry - ZNF415 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat brain using ZNF415 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - ZNF415 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using ZNF415 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - ZNF415 Polyclonal Antibody

ZNF415 Polyclonal Antibody



Catalog Number: 162337

Immunohistochemistry of paraffin-embedded

mouse kidney using ZNF415 antibody at dilution of 1:100 (40x lens).