

# Anti-beta Catenin Rabbit antibody



Catalog Number: 167337

## Product name

Rabbit IgG

Anti-beta Catenin Rabbit antibody

## Applications

## Specificity

WB, ICC/IF, FC

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/1000

beta Catenin Rabbit polyclonal antibody

ICC/IF: 1/50

## Preparation


FC: 1/50

Antigen: A synthetic peptide of human beta Catenin

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

Western blot detection of beta Catenin in C6 cell lysates using beta Catenin Rabbit pAb(1:1000 diluted).Predicted band size:86kDa.Observed band size:86kDa.

## Ig Type