

# NR2F6 antibody



Catalog Number: 107337

## Product name

NR2F6 antibody

## Specificity

Human; other species not tested.

## Antibody description

NR2F6 Mouse Monoclonal antibody. Positive IHC detected in human kidney tissue. Positive WB detected in MCF-7 cells, MCF7 cells. Observed molecular weight by Western-blot: 42 kDa

## Preparation

This antibody was obtained by immunization of NR2F6 recombinant protein (Accession Number: NM\_005234). Purification method: .

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human kidney using Catalog No:107337(NR2F6 antibody) at dilution of 1:50 (under 40x lens)



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107337(NR2F6 Antibody) at dilution of 1:1000