

# ATP11B antibody

Catalog Number: 108337

## Product name

ATP11B antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ATP11B Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Observed molecular weight by Western-blot: 12kd

## Preparation

This antibody was obtained by immunization of ATP11B recombinant protein (Accession Number: BC033880). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:108337(ATP11B antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:108337(ATP11B antibody) at dilution of 1:100 (under 40x lens)