

# PCGF6 antibody



Catalog Number: 113737

## Product name

PCGF6 antibody

## Specificity

Human; other species not tested.

## Antibody description

PCGF6 Rabbit Polyclonal antibody. Positive WB detected in PC-13 cells. Observed molecular weight by Western-blot: 39 kDa

## Preparation

This antibody was obtained by immunization of PCGF6 recombinant protein (Accession Number: NM\_001011663). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



PC-13 cells were subjected to SDS PAGE followed by western blot with Catalog No:113737(PCGF6 Antibody) at dilution of 1:1500