

# CTTNBP2NL antibody



Catalog Number: 109737

## Product name

CTTNBP2NL antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CTTNBP2NL Rabbit Polyclonal antibody. Positive IP detected in MCF-7 cells. Positive WB detected in MCF-7 cells, MDA-MB-453s cells, mouse skin tissue. Observed molecular weight by Western-blot: 70-75 kDa

## Preparation

This antibody was obtained by immunization of CTTNBP2NL recombinant protein (Accession Number: XM\_047425362). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:109737(CTTNBP2NL Antibody) at dilution of



IP Result of anti-CTTNBP2NL (IP:Catalog No:109737, 4ug; Detection:Catalog No:109737 1:600) with MCF-7 cells lysate 3200ug.