

POLR2L antibody



Catalog Number: 114047

Product name

POLR2L antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

POLR2L Rabbit Polyclonal antibody. Positive WB detected in Raji cells, HeLa cells. Observed molecular weight by Western-blot: 8 kDa

Preparation

This antibody was obtained by immunization of POLR2L recombinant protein (Accession Number: NM_021128). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:114047(POLR2L antibody) at dilution of 1:400