

CD1A antibody



Catalog Number: 107047

Product name

CD1A antibody

Specificity

Human; other species not tested.

Antibody description

CD1A Mouse Monoclonal antibody. Positive WB detected in Raji cells. Observed molecular weight by Western-blot: 49 kDa

Preparation

This antibody was obtained by immunization of CD1A recombinant protein (Accession Number: BC031645). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:107047(CD1A antibody) at dilution of 1:300