

CD133 antibody



Catalog Number: 109047

Product name

CD133 antibody

Immunogen

[Human CD133 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CD133 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, mouse pancreas tissue. Observed molecular weight by Western-blot: 115kd-120kd

Preparation

This antibody was obtained by immunization of CD133 recombinant protein (Accession Number: NM_001145847). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:109047(PROM1 antibody) at dilution of 1:500