

H11N2 HA (His Tag) recombinant protein



Catalog Number: 500447

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/thick-billed murre/Newfoundland/031/2007(H11N2)) hemagglutinin (ADD21478.1) (Met1-Lys528), termed as HA, was expressed with a polyhistidine tag at the C-terminus.

Organism

H11N2

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 17

Molecular Mass

The recombinant hemagglutinin of Influenza A virus (A/thick-billed

murre/Newfoundland/031/2007(H11N2)) consists 523 amino acids and predicts a molecular mass of 58.7 kDa.

Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, pH 7.4, 10 % glycerol.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H11N2 (A/thick-billed murre/Newfoundland/031/2007) Hemagglutinin / HA Protein (His Tag)