

TSPAN2 antibody

Catalog Number: 116447

Product name

TSPAN2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TSPAN2 Rabbit Polyclonal antibody. Positive IP detected in HL-60 cells. Positive WB detected in HL-60 cells. Observed molecular weight by Western-blot: 30-35 kDa

Preparation

This antibody was obtained by immunization of TSPAN2 recombinant protein (Accession Number: NM_005725). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:116447(TSPAN2 antibody) at dilution of 1:500



IP Result of anti-TSPAN2 (IP:Catalog No:116447, 4ug; Detection:Catalog No:116447 1:500) with HL-60 cells lysate 960ug.