

Human FAS/APT1/FAS1/TNFRSF6 (Fc tag) recombinant protein



Catalog Number: 518447

General Information

Gene Name Synonym

Tumor necrosis factor receptor superfamily member 6, Apo-1 antigen, Apoptosis-mediating surface antigen FAS, FASLG receptor, APT1, FAS1, TNFRSF6 and FAS

Protein Construction

A DNA sequence encoding the Human FAS/APT1/FAS1/TNFRSF6 Gln26-Asn173 is expressed with a Fc tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of PBS, pH 7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human FAS/APT1/FAS1/TNFRSF6 (Fc tag)
recombinant protein