

Anti-ALMS1 antibody

Catalog Number: 175847

Product name

Anti-ALMS1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to ALMS1

Preparation

This antigen of this antibody was synthesized peptide derived from alms1 . at aa range: 1530-1610

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.: Western Blot analysis of various cells using ALMS1 Polyclonal Antibody