

# AHSA1 antibody

Catalog Number: 107847

## Product name

AHSA1 antibody

WB: 1:200-1:2000

IP: 1:200-1:2000

## Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

IF: 1:20-1:200

## Antibody description

AHSA1 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas tissue. Positive IF detected in A431 cells. Positive IP detected in mouse brain tissue. Positive WB detected in A549 cells, A431 cells, HEK-293 cells, HeLa cells, mouse brain tissue. Observed molecular weight by Western-blot: 38 kDa

## Validations

### Preparation

This antibody was obtained by immunization of AHSA1 recombinant protein (Accession Number: NM\_012111). Purification method: Antigen affinity purified.

### Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

### Storage

Store at -20°C. DO NOT ALIQUOT

### Clonality

Polyclonal

### Ig Type


Rabbit IgG

### Applications


ELISA, IHC, IF, WB, IP

### Dilutions

Recommended Dilution:



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:107847(AHSA1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:107847(AHSA1 antibody) at dilution of 1:50 (under 40x lens)

