

Anti-RAB5C Rabbit antibody



Catalog Number: 168847

Product name

Anti-RAB5C Rabbit antibody

Specificity

Human

Antibody description

RAB5C Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human RAB5C

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

1:1000

1:20-1:100

1:50

Validations



Western blot detection of RAB5C in K562,C6,3T3,HeLa cell lysates using RAB5C Rabbit mAb(1:1000 diluted).Predicted band size:24kDa.Observed band size:24kDa.