

PPP1R11-Specific antibody



Catalog Number: 114147

Product name

PPP1R11-Specific antibody

Specificity

human, mouse; other species not tested.

Antibody description

PPP1R11-Specific Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A431 cells, mouse liver tissue. Observed molecular weight by Western-blot: 25-30kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_021959). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114147(PPP1R11-Specific antibody) at dilution of 1:100