

Human PSG2/PSBG2 (His tag) recombinant protein



Catalog Number: 519147

General Information

Gene Name Synonym

Pregnancy specific beta-1-glycoprotein 2, PSBG2, Pregnancy-specific beta-1 glycoprotein E, pregnancy-specific beta-1-glycoprotein 2, pregnancy-specific beta-1-glycoprotein 7, PSG2

Protein Construction

A DNA sequence encoding the Human PSG2/PSBG2 Gln35-Thr335 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human PSG2/PSBG2 (His tag) recombinant protein