

# SPATA1 antibody



Catalog Number: 115547

## Product name

SPATA1 antibody

## Immunogen

[Human SPATA1 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

SPATA1 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 40-50 kDa

## Preparation

This antibody was obtained by immunization of SPATA1 recombinant protein (Accession Number: BC127634). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115547(SPATA1 Antibody) at dilution of 1:300