

# UCP3 antibody

Catalog Number: 116547

## Product name

UCP3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

UCP3 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in mouse skeletal muscle tissue, SGC-7901 cells. Observed molecular weight by Western-blot: 30-34 kDa

## Preparation

This antibody was obtained by immunization of UCP3 recombinant protein (Accession Number: BC008392). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116547(UCP3 antibody) at dilution of 1:100



(under 10x lens)

mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:116547(UCP3 Antibody) at dilution of 1:1000