

# ZAP70 antibody



Catalog Number: 107547

---

## Product name

ZAP70 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ZAP70 Mouse Monoclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 70 kDa

## Preparation

This antibody was obtained by immunization of ZAP70 recombinant protein (Accession Number: BC039039). Purification method: .

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107547(ZAP70 antibody) at dilution of 1:1000