

RXRB antibody

Catalog Number: 114947

Product name

RXRB antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RXRB Rabbit Polyclonal antibody. Positive IP detected in MCF-7 cells. Positive WB detected in NIH/3T3 cells, MCF7 cells. Observed molecular weight by Western-blot: 57 kDa

Preparation

This antibody was obtained by immunization of RXRB recombinant protein (Accession Number: NM_021976). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:114947(RXRB antibody) at dilution of 1:500



IP Result of anti-RXRB (IP:Catalog No:114947, 4ug; Detection:Catalog No:114947 1:500) with MCF-7 cells lysate 1200ug.