

Anti-DDX54 antibody



Catalog Number: 175947

Product name

Anti-DDX54 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to DDX54

Preparation

This antigen of this antibody was synthesized peptide derived from ddx54 . at aa range: 570-650

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations



Fig1:; Western Blot analysis of various cells using DDX54 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit.