

# Anti-NF-kB p65 Rabbit antibody



Catalog Number: 166947

## Product name

Anti-NF-kB p65 Rabbit antibody

WB: 1/1000

ICC/IF: 1/20-1/100

## Specificity

Human

FC: 1/20

## Antibody description

NF-kB p65 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human NF-kB p65

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


WB, ICC/IF, FC

## Dilutions

## Validations



Western blot detection of NF-kB p65 in HeLa cell lysates using NF-kB p65 Rabbit pAb(1:1000 diluted). Predicted band size:60KDa. Observed



band size:65KDa.

Western blot detection of NF-kB p65 in Jurkat, C6 cell lysates using NF-kB p65 Rabbit pAb(1:1000 diluted). Predicted band size:60kDa. Observed band size:65kDa.