

# CBX5 antibody

Catalog Number: 108947

## Product name

CBX5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CBX5 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A431 cells, HEK-293 cells, human kidney tissue, K-562 cells, MCF7 cells, MCF-7 cells. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 25-30 kDa

## Preparation

This antibody was obtained by immunization of CBX5 recombinant protein (Accession Number: NM\_001127321). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:108947(CBX5 antibody) at dilution of 1:500



IP Result of anti-CBX5 (IP:Catalog No:108947, 4ug; Detection:Catalog No:108947 1:500) with HEK-293 cells lysate 1120ug.