

# Mouse IL-17F recombinant protein



Catalog Number: 501247

## General Information

### Protein Construction

A DNA sequence encoding the mouse Il17f (NP\_665855.2) (Arg21-Ala153) was expressed with an initial Met.

### Organism

Mouse

### Expression Host

E. coli

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE.

### Endotoxin

Please contact us for more information.

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

Met

### Molecular Mass

The recombinant mouse Il17f consists 134 amino

acids and predicts a molecular mass of 15 kDa.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Il17f protein SDS-PAGE