

EEA1 antibody



Catalog Number: 107247

Product name

EEA1 antibody

Specificity

Human; other species not tested.

Antibody description

EEA1 Mouse Monoclonal antibody. Positive IF detected in HeLa cells.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_003566). Purification method: Thiophilic affinity chromatograph.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgM

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:20-1:200

Validations



Immunofluorescent analysis of HeLa cells using (EEA1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L)