

LSM2 Polyclonal Antibody

Catalog Number: 165647

Product name

LSM2 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to LSM2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human LSM2 (NP_067000.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Western blot - LSM2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using LSM2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 30s.



Immunohistochemistry - LSM2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using LSM2 antibody at



dilution of 1:100 (40x lens).

LSM2 Polyclonal Antibody



Catalog Number: 165647

Immunofluorescence - LSM2 Polyclonal Antibody
Immunofluorescence analysis of U2OS cells using

LSM2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.