

CCAR2 Monoclonal Antibody



Catalog Number: 166647

Product name

CCAR2 Monoclonal Antibody

Specificity

Human, Mouse, Rat, Monkey

Antibody description

Monoclonal antibody to CCAR2

Preparation

Antigen: Recombinant protein of human CCAR2

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG

Applications

WB, ICC, IP


Dilutions

WB 1:200 - 1:1000

ICC 1:50 - 1:200

IP 1:20 - 1:50

Validations



Western blot - CCAR2 Monoclonal Antibody

Western blot analysis of extracts of various cell lines, using CCAR2 antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 20s.