

# Ubiquilin 1 antibody



Catalog Number: 116647

## Product name

Ubiquilin 1 antibody

## Specificity

Human, Rat; other species not tested.

## Antibody description

Ubiquilin 1 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue, L02 cells, rat heart tissue, rat liver tissue, Transfected HEK-293 cells. Positive IF detected in L02 cells. Observed molecular weight by Western-blot: 62-71 kDa

## Preparation

This antibody was obtained by immunization of Ubiquilin 1 recombinant protein (Accession Number: NM\_053067). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

## Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116647(Ubiquilin 1 Antibody) at dilution of 1:300



Immunofluorescent analysis of L02 cells using Catalog No:116647(Ubiquilin 1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)