

# Anti-TFAP2C/AP2 gamma Rabbit antibody



Catalog Number: 167647

## Product name

Anti-TFAP2C/AP2 gamma Rabbit antibody

## Specificity

Human

## Antibody description

TFAP2C/AP2 gamma Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human AP2 gamma/TFAP2C

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC

## Validations



Western blot detection of AP2 gamma/TFAP2C in HeLa,A549 using AP2 gamma/TFAP2C Rabbit mAb(1:1000 diluted)