

C12orf54 antibody

Catalog Number: 108647

Product name

C12orf54 antibody

Specificity

Human; other species not tested.

Antibody description

C12orf54 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human brain tissue.

Preparation

This antibody was obtained by immunization of C12orf54 recombinant protein (Accession Number: NM_152319). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:108647(C12orf54 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:108647(C12orf54 Antibody) at dilution of 1:50 (under 40x lens)