

Human AGRP (Fc Tag) recombinant protein



Catalog Number: 503347

General Information

Gene Name Synonym

Agouti-related protein

Protein Construction

A DNA sequence encoding the human AgRP (NP_001129.1) (Met 1-Thr 132) was expressed with the fused Fc region of human IgG1 at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 97 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ala 21

Molecular Mass

The recombinant human AgRP/Fc is a disulfide-linked homodimeric protein. The reduced monomer consists of 350 amino acids and predicts a molecular mass of 39 kDa also estimated by SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.2

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human AgRP Protein (Fc Tag) SDS-PAGE