

Phospho-GRIN1-S896 pAb



Catalog Number: 166347

Product name

Phospho-GRIN1-S896 pAb

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to Phospho-GRIN1-S896

Preparation

Antigen: A phospho specific peptide corresponding to residues surrounding S896 of human GRIN1

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - Phospho-GRIN1-S896 pAb

Western blot analysis of extracts from mouse Brain tissue, using phospho-GRIN1-S896 antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% BSA.