

Trappc9,NIBP antibody



Catalog Number: 116347

Product name

Trappc9,NIBP antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Trappc9,NIBP Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, MCF7 cells. Observed molecular weight by Western-blot: 130-139 kDa

Preparation

This antibody was obtained by immunization of Trappc9,NIBP recombinant protein (Accession Number: NM_180662). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:116347(Trappc9,NIBP antibody) at dilution of 1:1000