

Human NKG2A / NKG2 / CD159A / KLRC1 (His Tag) recombinant protein



Catalog Number: 501747

General Information

Gene Name Synonym

CD159 antigen-like family member A; NK cell receptor A; NKG2-A/B-activating NK receptor

Protein Construction

A DNA sequence encoding the mature form of human KLRC1 (NP_002250) (Arg100-Leu233) was expressed, with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 93 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human human KLRC1 consists of 154 amino acids and predicts a molecular mass of 17.8 KDa. It migrates as an approximately 30-40 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human NKG2A / CD159A / KLRC1 Protein (His Tag) SDS-PAGE