

PLOD3 Polyclonal Antibody



Catalog Number: 163747

Product name

PLOD3 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to PLOD3

Preparation

Antigen: A synthetic peptide of human PLOD3

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:200 - 1:2000

Validations



Western blot - PLOD3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PLOD3 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.