

Anti-ARID3B antibody



Catalog Number: 105747

Product name

Anti-ARID3B antibody

Specificity

Human ARID3B

Antibody description

Rabbit polyclonal to ARID3B

Preparation

Produced in rabbits immunized with A synthetic peptide corresponding to the N-terminus of the human ARID3B , and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 10-15 µg/mL

Validations



Lanes	A
Items	
Sample (whole cell lysate)	293T
Sample Volume (µg/lane)	30
Gel	13%SDS PAGE reducing gel
Recommended Concentration	10-15 µg/ml
Secondary Antibody	Dylight 800 labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.
Developed using Odyssey imaging system.	
Explanation	Predicted band size : 61 kDa Observed band size : 61 kDa Additional bands at : 40 kDa (We are unsure as to the identity of these extra bands. They are possible non-specific bindings).

ARID3B Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot