

# Anti-PBK/SPK Rabbit antibody



Catalog Number: 167057

## Product name

Rabbit IgG

Anti-PBK/SPK Rabbit antibody

## Applications

## Specificity

WB, IHC-P, IP

Human

## Dilutions

## Antibody description

WB: 1/5000-1/20000

PBK/SPK Rabbit polyclonal antibody

IHC: 1/20-1/100

## Preparation

IP: 1/20

Antigen: A synthetic peptide of human PBK

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage


Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

## Validations



Western blot detection of PBK in K562, HeLa cell lysates using PBK Rabbit pAb(1:1000 diluted). Predicted band size: 36kDa. Observed band size: 40kDa.