

Anti-His-tag (1H3) Mouse antibody



Catalog Number: 169057

Product name

WB

Anti-His-tag (1H3) Mouse antibody

Dilutions

Antibody description

WB:1/500-1/2000

His-tag (1H3) Mouse monoclonal antibody

Validations

Preparation

Antigen: Synthetic peptide corresponding to fragment (HHHHHH) His Tag

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications



Western blot detection of His-tag in various bacterial lysates, which containing His-tagged proteins (MW: 30 kDa, 45 kDa, 72 kDa, 44 kDa and 40 kDa respectively) using anti-his antibody (1:1000 diluted).