

## General Information

### Gene Name Synonym

IL17Rhom; Interleukin-17 receptor-like protein;  
Sef homolog

### Protein Construction

A DNA sequence encoding the extracellular domain of human IL17RD (NP\_060033.3) precursor (Met 1-Arg 299) was expressed with a C-terminal polyhistidine tag.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 97 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Cys 17

### Molecular Mass

The secreted recombinant human IL17RD comprises 294 amino acids with a predicted molecular mass of 33.5 kDa. As a result of glycosylation, the apparent molecular mass of rhIL17RD is approximately 55-60 kDa in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4  
1. 5 % trehalose and mannitol are added as protectants before lyophilization.  
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human IL17RD Protein (His Tag) SDS-PAGE