

# Anti-ADAMTS18 antibody



Catalog Number: 176457

## Product name

Anti-ADAMTS18 antibody

## Specificity

Human, Mouse

## Antibody description

Rabbit polyclonal antibody to ADAMTS18

## Preparation

This antigen of this antibody was synthesized peptide derived from adamts-18 . at aa range: 1030-1110

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P

## Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

## Validations




Fig1:: Western blot analysis of 3T3 using ADAMTS-18 antibody.. Secondary antibody [catalog#: HA1001) was diluted at 1:20000