

# Anti-COX2/cyclooxygenase2 Rabbit antibody



Catalog Number: 167457

## Product name

Anti-COX2/cyclooxygenase2 Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

WB, IHC-P, ICC/IF, ELISA

## Antibody description

COX2/cyclooxygenase2 Rabbit polyclonal antibody

## Dilutions

WB: 1/500 - 1/2000

IHC: 1/100 - 1/300

IF: 1/200 - 1/1000

## Preparation

Antigen: The antiserum was produced against synthesized peptide derived from human Cox2. AA range:555-604

ELISA: 1/20000

## Validations

## Formulation


Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal



Western Blot analysis of COx-2 in A549 cells using Cox-2 Polyclonal Antibody at dilution 1/2000