

# Anti-Rad51 (4H10) Mouse antibody



Catalog Number: 168457

## Product name

Anti-Rad51 (4H10) Mouse antibody

## Specificity

Human, Mouse, Rat, Monkey

## Antibody description

Rad51 (4H10) Mouse monoclonal antibody

## Preparation

Antigen: Recombinant protein of human Rad51.

## Formulation

0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of Rad51 in CHO-K1, COS7, C6, 3T3 and HeLa cell lysates using Rad51 mouse mAb (1:1000 diluted). Predicted band size: 37KDa. Observed band size: 37KDa.