

BMF antibody



Catalog Number: 108457

Product name

BMF antibody

Specificity

Human; other species not tested.

Antibody description

BMF Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, HEK-293 cells. Observed molecular weight by Western-blot: 24kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001398507). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:117192(BMF antibody) at dilution of 1:300