

Anti-AQP12A antibody



Catalog Number: 175857

Product name

Anti-AQP12A antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to AQP12A

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human aqp12. aa range:231-280

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations



Fig1.; Western Blot analysis of various cells using AQP12 Polyclonal Antibody diluted at 1: 1000