

Human CGB5/CG-beta (Fc Tag) recombinant protein



Catalog Number: 503157

General Information

Gene Name Synonym

Chorionic gonadotrophin chain beta

Protein Construction

A DNA sequence encoding the human CGB5 (NP_000728.1) (Met1-Gln165) was expressed with the Fc region of mouse IgG1 at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Ser 21

Molecular Mass

The recombinant human CGB5/mFc comprises 379 amino acids and has a predicted molecular mass of 41.9 kDa. The apparent molecular mass of the monomer is approximately 56 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human CGB5 Protein (Fc Tag) SDS-PAGE