

TNPO2 antibody



Catalog Number: 116157

Product name

TNPO2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TNPO2 Rabbit Polyclonal antibody. Positive IF detected in MCF-7 cells. Positive WB detected in HeLa cells, HEK-293 cells, human brain tissue, Jurkat cells, K-562 cells. Observed molecular weight by Western-blot: 100kd

Preparation

This antibody was obtained by immunization of TNPO2 recombinant protein (Accession Number: NM_001382236). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116157(TNPO2 antibody) at dilution of 1:500



Immunofluorescent analysis of MCF-7 cells, using TNPO2 antibody Catalog No:116157 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).