

# ZNF3 antibody



Catalog Number: 117157

## Product name

ZNF3 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

ZNF3 Rabbit Polyclonal antibody. Positive WB detected in mouse ovary tissue, A2780 cells, mouse testis tissue. Observed molecular weight by Western-blot: 47-51 kDa

## Preparation

This antibody was obtained by immunization of ZNF3 recombinant protein (Accession Number: NM\_032924). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse ovary tissue were subjected to SDS PAGE followed by western blot with Catalog No:117157(ZNF3 antibody) at dilution of 1:1000